# promed EMT-6

## Model: LT1021A

Digitales, elektrisches Schmerztherapie- und Muskelstimulations-Gerät (TENS & EMS)  
(TENS: Transkutane Elektrische Nerven-Stimulation) (EMS: Elektrische Muskel-Stimulation)

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**LCD Display Parts:**

- 1.** Intensity of channel
- 2.** Treatment program
- 3.** Display of therapeutic mode
- 4.** Display of pulse rate
- 5.** EMS waveform of working time
- 6.** Display of pulse width
- 7.** EMS waveform of rest time
- 8.** EMS waveform of ramp up and ramp down time
- 9.** Time symbol
- 10.** Parameters and treatment time
- 11.** Low battery indicator
- 12.** Output waveform

|  |          |   |            |
|--|----------|---|------------|
| GETTING STARTED                              | p. 21    | NORMALIZED SYMBOLS                        | S. 201/202 |
| CONTENTS                                     | p. 22    | CE DIRECTIVES                             | S. 203/204 |
| SYSTEM COMPONENTS                            | p. 23    | TECHNICAL SPECIFICATIONS                  | S. 205/206 |
| DESCRIPTION OF THE UNIT                      | p. 24    | TRANSPORT/STORAGE/OPERATING<br>CONDITIONS | S. 207/208 |
| INTRODUCTION                                 | p. 25/26 | ACCESSORIES                               | S. 229/230 |
| INDICATIONS / CONTRAINDICATIONS              | p. 27/28 | CLINICAL RECOMANDATIONS                   | S. 231/232 |
| WARNINGS / PRECAUTIONARY MEASURES            | p. 29/30 | ANNEX A PLACEMENT OF ELECTRODES           | S. 233     |
| OPERATING THE UNIT                           | p. 31-34 | GUARANTEE CARD                            | S. 234     |
| PROGRAM                                      | p. 35/36 | MANUFACTURER + FURTHER INFORMATION        | S. 0       |
| SAFETY INSTRUCTIONS                          | p. 37    |   |            |
| TROUBLESHOOTING                              | p. 38    |   |            |
| CARE, MAINTENANCE AND STORAGE                | p. 39    |   |            |
| DISPOSAL                                     | p. 39    |   |            |
| GUARANTEE CONDITIONS<br>AND CUSTOMER SERVICE | p. 40    |   |            |

DE

EN

FR

IT

ES

NL

RU


PL

SE

FI

Your **Promed EMT-6 - package** contains the following equipment:

- |                             |       |
|-----------------------------|-------|
| 1. EMT-6 Unit               | 1 PCS |
| 2. Cable                    | 2 PCS |
| 3. Electrode pads (40*40mm) | 4 PCS |
| 4. User manual              | 1 PCS |
| 5. AAA battery              | 4 PCS |
| 6. Carrying case            | 1 PCS |

We reserve the right to make changes for the purpose of technical improvement without prior notice.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

Please read the user manual thoroughly before using the device for the first time.

Thank you for purchasing the **Promed EMT-6**. You have purchased a high-quality product designed for personal hygiene and wellbeing purposes. **Promed** is a leading company with decades of experience in personal hygiene, wellness and health areas.

The **Promed EMT-6** unit has been designed and manufactured in accordance with the Medical Device Directives 93/42/EEC and 2007/47/ EC to guarantee quality during use. The device may be used once this user manual has been read.

We, as the manufacturer, cannot be made liable in any way for injury or damages to people or objects that arise from failure to comply with this user manual. We wish you a lot of enjoyment with your new **Promed EMT-6**. We would like to familiarize you with the therapy unit in the following sections. Please read the user manual thoroughly before using the device for the first time.

You can receive professional advice wherever **Promed** products are sold or get in touch with us if you have any questions. We can provide you with the name of the representative responsible for you.

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Further informations can be found on our homepage:  
**[www.promed.de](http://www.promed.de)**

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

**Promed EMT-6** electrical stimulator is a portable electrotherapy device featuring two therapeutic modes: Transcutaneous Electrical Nerve Stimulation (TENS), and Electrical Muscle Stimulation (EMS), which are used for pain relief and electrical muscle stimulation. The stimulator sends gentle electrical current to underlying nerves and muscle groups via electrodes applied on the skin. The parameters of the device are controlled by the press buttons. The intensity level is adjustable according to the needs of patients. Before using, please read all the instructions in this user manual carefully and keep it safe for future reference.

### What is pain?

Pain is the body's own early warning system. Feeling pain is important, as it indicates an abnormal condition within the body and serves as a warning to us before any further damage or injury occurs. However, long-term, persistent pains – often known as chronic pain – serve no obvious useful purpose once a diagnosis has been made. TENS was developed to soothe or eliminate certain types of chronic and acute pain.

### We differentiate between two types of pain:

#### • Acute pain

as a chief symptom can often help the doctor with diagnosis and the acute pain has a protective function for the patient.

#### • Chronic pain

can often become part of the illness in itself. A patient suffering from chronic pain will often suffer for years and experiences changes in his/her personality structure.

### EXPLANATION OF TENS

Transcutaneous Electrical Nerve Stimulation (TENS) is a non-invasive, drug free method of controlling pain. TENS uses tiny electrical impulses sent through the skin to nerves to modify your pain perception. TENS does not cure any physiological problem; it only helps control the pain. TENS does not work for everyone; however, in most patients it is effective in reducing or eliminating the pain, allowing for a return to normal activity.

### How does TENS work?

**Promed EMT-6** works by passing harmless electrical signals into the body from the electrode pads. This relieves pain in two ways:

- Firstly, it blocks the body's pain signals which are normally transmitted from the area of damage through the nerve fibres to the brain - TENS interrupts these pain signals.
- Secondly, TENS stimulates the body's production of endorphins – its own natural painkillers.

### Possible TENS applications

Generally the **Promed EMT-6** shall be used for the following medical indications or the following complaints as a TENS unit for treatment:

- For symptomatic relief of chronic intractable pain
- For post traumatic pain (acute pain occurring)
- For post surgical pain (pain caused by surgery)

## EXPLANATION OF EMS

Electrical Muscle Stimulation (EMS) is an internationally accepted and proven way of treating muscular injuries. It works by sending electronic pulses to the muscle needing treatment; this causes the muscle to exercise passively. It is a product derived from the square waveform, originally invented by John Faraday in 1831. Through the square wave pattern it is able to work directly on muscle motor neurons. This device has low frequency and in conjunction with the square wave pattern allows direct work on muscle groupings. This is being widely used in hospitals and sports clinics for the treatment of muscular injuries and for the re-education of paralyzed muscles, to prevent atrophy in affected muscles and improving muscle tone and blood circulation.

### How does EMS work?

The EMS units send comfortable impulses through the skin that stimulate the nerves in the treatment area. When the muscle receives this signal it contracts as if the brain has sent the signal itself. As the signal strength increases, the muscle flexes as in physical exercise. Then when the pulse ceases, the muscle relaxes and the cycle is repeated.

The goal of electrical muscle stimulation is to achieve contractions or vibrations in the muscles. Normal muscular activity is controlled by the central and peripheral nervous systems, which transmit electrical signals to the muscles. EMS works similarly but uses an external source (the stimulator) with electrodes attached to the skin for

transmitting electrical impulses into the body. The impulses stimulate the nerves to send signals to a specifically targeted muscle, which reacts by contracting, just as it does with normal muscular activity.

### Possible EMS application

Generally the **Promed EMT-6** shall be used for the following medical indications or the following complaints as a EMS unit for treatment:

- For relaxation of muscle spasm
- For increase of blood circulation
- For prevention of retardation of disuse atrophy
- For muscle re-education

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

*Be sure to consult a doctor before using the unit!*

*In the USA, the sale of this device is restricted by law to doctors, or subject to a doctor's prescription.*

Please carefully read and understand the following warnings and cautions to ensure the safe and correct use of this device and to prevent injury.

### Indications

The device is designed to be used for temporary relief of pain, including the acute and chronic pain relief.

### Contraindications

Do not use this device without consulting a doctor if you have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device. Such use could cause electric shock, burns, electrical interference, or death.

If you apply to one of the following cases always be sure to consult a doctor before you apply the **Promed EMT-6** and clarify the use of the device with him.

- In case of persistent pain symptoms, despite therapy
- Taking a strong analgesic or local anesthetic
- In Infectious Diseases
- By circulatory disorder (thrombosis and embolism)
- For sensory disturbances (numbness)
- When used in infants, toddlers and children
- In pregnancy
- With Psychosis
- If bleeding tendency
- With Cancer
- At extreme sensitivity to or fear of electricity
- In patients with metallic implants
- For heart problems, especially cardiac arrhythmias
- Prior to each electrode placement, that apply current to the carotid sinus region (anterior neck)
- Prior to each electrode placement, that apply current transcebrally (through the head)
- If undiagnosed pain symptoms
- For treatment of the eyelids

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

- Serious arterial circulatory problems (embolism) in the lower limbs
- The device should not be used for symptomatic local pain relief unless etiology is established or unless a pain syndrome has been diagnosed
- When cancerous lesions are present in the treatment area
- At swollen, infected, inflamed areas or skin eruption (e.g. phlebitis, thrombophlebitis, varicose veins, etc.)
- On-existence of a demand-type cardiac pacemaker or any implanted defibrillator
- In areas of the body with poorly enervated areas
- With epilepsy
- For abdominal or inguinal hernia
- Do not apply stimulation over your neck because this could cause severe muscle spasms resulting in closure of your airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure;

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

**⚠ Warnings**

- TENS devices must be kept out of the reach of children. Avoid inhalation or swallowing of small parts. And the cable may cause strangulation.
- The safety of using TENS devices during pregnancy or delivery has not been proven.
- TENS is not effective in treating pain of the central nervous system (headaches).
- Should TENS treatment fail to work or cause discomfort, stimulation should be discontinued until a doctor has reassessed the situation.
- Always switch the TENS unit off before attaching or removing the electrodes.
- Never position the electrodes over the eyes, in the mouth or within the body.
- TENS devices do not possess any healing properties.
- TENS is a method for treating symptoms and as such only suppresses the perception of pains that would otherwise serve as a defence and warning mechanism for the body.
- Do not apply stimulation while driving, operating machinery, or during any activity in which electrical stimulation can put you at risk of injury.
- If you are in the care of a physician, consult with your physician before using this device.
- If you have had medical or physical treatment for your pain, consult with your physician before using this device.
- If your pain does not improve, becomes more than mild, or continues for more than five days, stop using the device and consult with your physician.
- Do not apply stimulation over your neck because this could cause severe muscle spasms resulting in closure of your airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Do not apply stimulation across your chest because the introduction of electrical current into the chest may cause rhythm disturbances to your heart, which could be lethal.
- Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use.
- Do not apply stimulation when in the bath or shower.
- Do not apply stimulation while sleeping.
- Do not use the device on children, if it has not been evaluated for pediatric use.
- Consult with your physician before using this device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals.
- Apply stimulation only to normal, intact, clean, healthy skin.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

### Precautionary measures / Negative reactions

#### Precautions

- TENS is not a substitute for pain medications and other pain management therapies.
- TENS devices have no curative value.
- TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism.
- The success of the therapy depends to a considerable extent on the patient's choice of physician and whether they are qualified for the treatment of patients suffering from pain.
- The long-term effects of electrical stimulation are unknown.
- You may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium.
- Use caution if you have a tendency to bleed internally, such as following an injury or fracture.
- Consult with your physician prior to using the device after a recent surgical procedure, because stimulation may disrupt the healing process.
- Use caution if stimulation is applied over the menstruating or pregnant uterus.
- Use caution if stimulation is applied over areas of skin that lack normal sensation.
- Use this device only with the leads, electrodes, and accessories recommended by the manufacturer.
- Do not modify the device or the electrodes without authorization of the manufacture. This could cause improper functioning.

- Because the effects of stimulation on the brain are unknown, do not use the device on opposite sides or your head.
- This device is designed for use by and on a single adult person.
- For hygienic reasons electrodes should not be shared.
- Do not modify electrodes (e.g. by cutting them). This increases the current density, which is potentially hazardous and effective current density beyond 2mA/cm<sup>2</sup> requires increased awareness.

#### Negative Reactions

- Long-term use can occasionally lead to skin irritation in the area where the electrodes were placed.
- The effectiveness of the treatment depends greatly upon the patient being treated by someone who is qualified in dealing with patients suffering from pain.
- Potential negative reactions include skin irritation and burns caused by the electrodes.
- You may experience skin irritation and burns beneath the stimulation electrodes applied to your skin.
- You may experience headache and other painful sensations during or following the application of electrical stimulation near your eyes and to your head and face.
- You should stop using the device and should consult with your physician if you experience adverse reactions from the device.

DE

EN

FR

IT

ES

NL

RU


PL

SE

FI

**Before you start:**

- Make sure the batteries are installed correctly.
- Connect the cable to the electrode pads and then connect to the unit.
- Attach the electrode pads to the area of pain.
- Do not use the electrode pads if it is scratched or damaged in any way.


**Check/Replace the battery**

1. Remove the Belt Clip (J) by sliding it down the unit.
2. Open the battery compartment cover (I).
3. Insert 4 batteries (type: AAA) into the battery compartment. Make sure you are installing the batteries correctly. Be sure to match the positive and negative ends of the batteries to the markings in the battery compartment of the device.
4. Recover the battery compartment (I) cover if you want to use it.
5. Bring on the Belt Clip (J) back by sliding it upwards the unit.

**Caution:**

1. Remove the batteries if the device is not in use for long periods of time.
2. Do not mix old and new batteries or different types of batteries.

3. Warning: If batteries leak and come into contact with the skin or eyes, wash immediately with copious amounts of water.
4. Batteries must be handled by an adult. Keep batteries out of the reach of children.
5. Only batteries of the same or equivalent type are recommended.
6. Do not use rechargeable batteries.
7. Remove exhausted batteries from the unit.
8. Dispose of batteries safely according to battery manufacturer's instructions.

**Connect the electrode pads to the cables:**

Hold the cable plug and insert into electrode pad connections. Make sure there are no bare metal of the pins exposed.

**Caution:**

Always use the cable that supplied with manufacturer or distributor and use the electrode pads with CE mark, or are legally marketed in the US under 510 (k) procedure.

**Connect cable to unit**

Before proceeding to this step, ensure the device is completely switched off. Hold the cable plug and insert into the socket (B) on the top of the unit.

**Caution:**

Do not insert the plug of the patient lead wire into any AC power supply socket.

**Place electrode pads on skin**

Apply electrode pads to the area where the pain is felt (refer to Annex A: placement of electrodes). Before applying electrodes, be sure the skin surface over which electrode pads are placed is thoroughly clean and dry. Make sure the electrode pads are pressed firmly to the skin and make good contact between the skin and the electrode pads. Place the electrode pads over the skin, attach them properly, firmly, and evenly.


**Caution:**

1. Before applying the electrode pads to the body, ensure the skin surface is clean and free from lotions or moisturizers.
2. Do not switch the device on when self-adhesive electrode pads are not positioned on the body.
3. Never remove the pads from the skin whilst the device is still switched on.
4. For hygiene reasons we recommend that the electrode pads are replaced every 30 days.
5. Please make sure the device is completely turned off before connecting.
6. It is recommended that, at minimum, 40mm x 40mm self-adhering based, square electrodes are used at the treatment area.
7. For hygienic reasons, each patient should use their own set of electrodes.

**Turn the device on:**

Push up the “Power ON/OFF” switch (A) on the right side of unit to turn the device on, the following screen will appear:


**Select the program:**

There are 11 (P03~P13) preset programs and 4 (P1, P2, S, A) user programs available for you to select. The details of the programs, please refer to section programs. Press the “MODE” button (E), you can select the therapeutic program you need. And the LCD will display like following:


### Set the parameter

TENS mode, in user programs (P1, P2), the user can set the pulse rate, pulse width and treatment time by press ["SET" button (F). The user LCD will be display the current parameter which you want to set like following:

EMS mode, in user programs (S, A), the user can set the pulse rate, pulse width, treatment time, ramp time, working time and rest time by press "SET" button (F). With programme 'A' (alternating) the programme initially proceeds on channel 1 before changing to channel 2. This changeover takes place during the full duration of the therapy. With programme 'S' (simultaneous), no changeover takes place. The programme proceeds on both channels simultaneously during the full duration of the therapy. The LCD will be display the current parameter which you want to set like following:

And then press "▼" (G) or "▲" (H) button to set the parameter. After the parameter is set, press the "SET" button (F) again to set next one.

**Remark:** The set range of pulse rate, pulse width, treatment time, ramp time, working time and rest time in user program, please refer to section programs (page 13-14).

### Restore to factory preset:

If you want to restore the user programs to factory reset, you can achieve by press the "▼" (G) button and turn the device on at the off state mode.

### Adjust the intensity and Start to treatment:


There are two intensity adjustment knobs (C) for each channel on the side of the unit. By turning the controls up, the output intensity will be increased, when the output intensity above 1 mA, the device will be started to work. The maximum output intensity is 100 mA.

### Remark:

The **Promed EMT-6** has special backlight function:

- In the normal state the backlight is green;
- If the output intensity is over 65 mA, the backlight will be changed blue and the user should be taken care.
- If the electrode pad connection in wrong way or not connect the electrodes, when the output intensity over 10 mA, the backlight will change to red and flash, and the intensity will be reseted automatically.

## Lock function

In order to use the device safely, after you adjusted the output intensity, please press the “SET” button (F) and “▲” button (H) at the same time, the lock function mark “🔑” (9) will show up on the LCD (The device will be locked). This is a safety feature to prevent accidental

changes to your settings and to prevent accidental increases the intensity level. If you want to unlock, Press the “SET” button (F) and “▲” button (H) at the same time, the lock function mark “🔑” (9) will disappear.

## Check the recording memory

1. The **Promed EMT-6** device can record 30 times treatment data. In the waiting or setting state, to check the memory data, press the “SET” button (F) and hold on for 5 seconds, the device will enter the memory state. First, the last therapeutic program recoding will display on the LCD, you can press the “SET” button (F) again to check the treatment parameter of this program (treatment time, pulse rate and pulse width).
2. Press “▼” (G) or “▲” (H) button to check other therapeutic program memory.
3. If you want to back the waiting or setting state, press the “MODE” button (E) or wait for 30 seconds without any operation in the panel.

4. To clear the memory, press the “SET” button (F) and hold on for 5 seconds, the symbol of “D” will be flash on the LCD to remind the user to delete the memory. You can press the “SET” button (F) again to delete the memory.

## Low battery indicator:

When the low power indicator (8) flashes, the batteries should be replaced with new ones as soon as possible. However, the device will continue to operate for several more hours. When the device power off, the current state will be saved automatically.

## Turn the device off

Push “Power ON/OFF” switch (A) down to turn off the device.

## Caution

If there is no operation in the panel in waiting state, you will hear a “DI” sound after one minute and hear two “DI” sounds after two minutes; and then the device will automatically turn into the power saving state after three minutes, LCD will close the display and you will hear three times “DI” sounds. You can press any button in the panel to activate the display.

There are 11 (P03~P13) preset programs and 4 (P1, P2, S, A) user programs available for you to select from. The Treatment time, pulse rate and pulse width can be set in TENS user program P1 and P2, and the pulse rate, pulse width, treatment time, ramp time, working time and rest time can be set in EMS user program S and A. The setting ranges of the parameter please refer to below. And the electrode pad positions please see Annex A Placement of electrodes.

### Programs for TENS

| Program                | Frequency      | Pulse width      | Waveform   | Treatment time               | Indication for use /<br>Electrode placement no.   |
|------------------------|----------------|------------------|------------|------------------------------|---|
| <b>user programs</b>   |                |                  |            |                              |   |
| [P1]                   | 20 - 110 Hz    | 50 - 200 $\mu$ s | Continuous | 1 min - 60 min<br>/ Continue | Neck pain 2<br>Shoulder pain 3<br>Elbow pain 7<br>Rheumatic pain 9<br>Lumbago 11,12<br>Menstrual pain 13<br>Phantom limb pain 14<br>Hip pain 16<br>Osteoarthritis pain<br>in the knee 18,19<br>wound healing 20 |
| [P2]                   | 10 Hz - 110 Hz | 200-100 $\mu$ s  | Modulation | 1 min - 60 min<br>/ Continue | Shoulder pain 3<br>Trapezius pain 10, 21<br>Lumbago 12<br>Thigh pain 22   |
| <b>preset programs</b> |                |                  |            |                              |   |
| P03                    | 110 Hz         | 50 $\mu$ s       | Continuous | 30 min                       | Facial pain 1<br>Neck pain 2  |
| P04                    | 4 Hz           | 200 $\mu$ s      | Continuous | 30 min                       | Postoperative or Chemotherapy-<br>induced Nausea 8  |

|     |  |             |                           |        |  |                      |
|-----|--|-------------|---------------------------|--------|--|----------------------|
| P05 | Burst rate: 2 Hz<br>Frequency fixed:<br>100 Hz | 200 µs      | Burst                     | 30 min | Cervical rhizopathy<br>Central pain<br>Sciatica<br>Knee pain | 4<br>5,6<br>15<br>19 |
| P06 | 110 Hz   | 200 µs      | Continuous                | 30 min | See program P1   |                      |
| P07 | 110 Hz   | 50 - 250 µs | Pulse width<br>Modulation | 30 min | Trapezius pain   | 10, 21               |
| P08 | 20-110 Hz                                      | 200 µs      | Pulse rate<br>modulation  | 30 min | Lumbago  | 11, 12               |

**Programs for EMS**

| Program                | Frequency   | Pulse width | Waveform   | Treatment time               | Work sec. | Rest sec. | Ramp sec. | Indication for use   |
|------------------------|-------------|-------------|------------|------------------------------|-----------|-----------|-----------|--|
| <b>preset programs</b> |             |             |            |                              |           |           |           |  |
| P09                    | 10 Hz       | 250 µs      | Continuous | 30 min                       | 3         | 6         | 2         | Urge and stress incontinence (electrode placement number 23/24)<br><br>Relaxation of muscle spasms;<br><br>Increase of blood flow circulation;<br><br>Prevention of disuse atrophy;<br><br>Muscle re-education |
| P10                    | 50 Hz       | 300 µs      | Continuous | 30 min                       | 5         | 10        | 1         |  |
| P11                    | 50 Hz       | 300 µs      | Continuous | 30 min                       | 5         | 15        | 1         |  |
| P12                    | 75 Hz       | 300 µs      | Continuous | 30 min                       | 5         | 10        | 1         |  |
| P13                    | 75 Hz       | 300 µs      | Continuous | 30 min                       | 5         | 15        | 1         |  |
| <b>user programs</b>   |             |             |            |                              |           |           |           |  |
| S                      | 10 - 110 Hz | 50 - 330 µs | Continuous | 1 min - 60 min<br>/ Continue | 1 - 30 s  | 1 - 60 s  | 1 - 6 s   |  |
| A                      | 10 - 110 Hz | 50 - 330 µs | Continuous | 1 min - 60 min<br>/ Continue | 1 - 30 s  | 1 - 60 s  | 1 - 6 s   |  |

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

**Pay attention to the following before first time operation:**

- The Promed EMT-6 device is only intended for domestic use.
- Scrupulous supervision is required when the device is being used in the vicinity of children.
- Never place or use the unit in a wet or moist environment.
-  Do not use the unit under water, e.g., under the shower.
- Never place the unit in direct contact with fire, gas or oxygen, or hot objects such as hot plates.
- Take every precaution possible to ensure that the unit does not fall down or get damaged in any other way.
- If problems do arise with the unit, take it immediately to get it repaired.
- Do not lubricate or wash the unit.

**⚠ Danger!**

- Do not place the unit in direct contact with water or other fluids.
- It is not permitted to use the unit outdoors.
- Never touch the unit with wet hands.
- Do not store the unit near a wash basin or bath tub as there is the danger that the unit may fall or be pulled into the wash basin or the bath tub.

**⚠ Warning!**

- Never leave the unit unattended when children or inexperienced persons are in the vicinity of the unit.
- Make sure that children do not play with the device.
- Only use the unit for the applications described in this manual.
- This unit is not designed for use by people (including children) who have limited physical, sensory, or intellectual abilities or who do not have sufficient experience and/or sufficient knowledge. This does not apply if these people are supervised by a person responsible for their safety or receive instructions from them on how to use the unit.
- Do not use the unit underneath blankets or cushions as excessive heat may cause fire, injury or an electrical shock.
- Do not carry the unit by its cable and do not use the cable as a handle.
- It is not permitted to use the unit outdoors.
- When you have finished treatment, place the intensity controllers **(2)** in the "OFF" position.
- Never let children play with the packaging material; there is a risk of suffocation.

| Problem   | Possible Cause   | Solution   |
|---|--|--|
| Displays fail to light up                       | Battery contact failure                                    | 1. Try fresh batteries.  |
|   |  | 2. Ensure batteries are inserted correctly. Check the following contacts: <ul style="list-style-type: none"> <li>• All contacts are in place.</li> <li>• All contacts are not broken.</li> </ul> |
| Stimulation weak or cannot feel any stimulation | Electrodes<br>1. Dried out or contaminated<br>2. Placement | Replace and re-connect   |
|   | Lead wires<br>Old/worn /damaged                            | Replace  |
|   | The intensity is too weak.                                 | Use a higher intensity level.  |
| Stimulation is uncomfortable.                   | Intensity is too high                                      | Decrease intensity.  |
|   | Electrodes are too close together                          | Reposition the electrodes.   |
|   | Damaged or worn electrodes or lead wires                   | Replace  |
|   | Electrode active area size is too small.                   | Replace electrodes with ones that have an active area no less than 16.0 cm <sup>2</sup> (4cm * 4cm).   |
|   | Mayn't operate the device according to the manual.         | Please check the manual before use.  |
| Intermittent output                             | Lead wires   | 1. Verify connection is secure. Insure firmly.   |
|   |  | 2. Turn down the intensity. Rotate lead wires in socket 90°. If still intermittent, replace lead wire.   |
|   |  | 3. If still intermittent after replacing lead wire, a component may have failed. Call the repair department.   |
|   |  | 4. Some programs will seem intermittent. This is expected.   |

|  |   |  |
|--|---|--|
| Stimulation is ineffective.                          | Improper electrode and applicator placement<br>Unknown. | Reposition electrode and applicator-Contact clinician.                                       |
| The skin becomes red and/or you feel a stabbing pain | Use the electrodes on the same site every time.         | Re-position the electrodes. If at any time you feel pain or discomfort stop use immediately. |
|  | The electrodes aren't stuck onto the skin properly.     | Ensure the electrode is stuck securely on the skin.  |
|  | The electrodes are dirty.                               | Clean the electrode pads with a damp, lint free cloth or replace new electrode pads.         |
| Output current stops during therapy                  | The surface of the electrode was scratched.             | Replace new electrode.   |
|  | The electrode pads come off the skin.                   | Turn off the device and stick the electrode pad firmly to the skin.                          |
|  | The cable is disconnected.                              | Turn off the device and connect the cable  |
|  | The power of the batteries has been exhausted.          | Please replace them with new batteries.  |

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

- Do not place the unit in direct contact with water or other fluids.
- Switch the unit off when it is not in use.
- Before putting the **Promed EMT-6** into storage for a longer period of time, remove the batteries from the compartment. Leaking batteries can damage the unit.
- Keep the unit and its accessories in a cool, dry place in the case provided.
- Handle the electrodes by their edges when removing them. To avoid damaging them, never pull at the electrode leads themselves.
- Do not make any sharp kinks in the connecting leads or electrodes.
- After use, stick the electrode pad onto the protective plastic film.
- Do not expose the device to direct sunlight and protect it against dirt and moisture.
- Never place any heavy objects on the device.
- You can clean your **Promed EMT-6** by wiping it over carefully with a cloth moistened with mild soapy water. You may also use isopropyl alcohol or a soap solution. Household detergents and cleaning products are not suitable.
- The unit must not be returned to your dealer or the manufacturer for inspection and recalibration. If you need such inspections or recalibrations because of your internal Quality Management System don't hesitate to contact Promed.
- If you should experience other problems, consult your dealer, returning the unit if necessary. Never attempt to repair a fault yourself.

Power tools, accessories, and packaging should be recycled in an environmentally suitable manner.

**Only for EU countries:** Do not throw power tools in with the domestic waste! In accordance with the European Directive 2012/19/EU on waste electrical and electronic equipment and its implementation into national law, nonfunctioning  power tools must be collected separately and recycled in an environmentally suitable manner.

This symbol indicates that disposal of this product in the domestic waste is not permitted within the EU. Waste devices contain valuable recyclable material that should be recycled. Recycling also prevents uncontrolled waste disposal from damaging the environment and human health. Therefore, please dispose of all waste devices using the appropriate collection systems or send the device to the place of purchase for disposal. The place of purchase will then recycle the device.

#### **Rechargeable batteries/batteries:**

Do not throw your rechargeable batteries/batteries in the domestic waste, in fire, or in water. Rechargeable batteries/batteries should be collected, recycled, or disposed of in an environmentally friendly manner.

#### **Only for EU countries:**

According to the directive 91/157/EEC, defect or used rechargeable batteries/batteries must be recycled. Waste rechargeable batteries/batteries can be directly handed to:

**Promed GmbH**, Lindenweg 11, D-82490 Farchant.

This product was produced with the utmost care and was thoroughly tested before it left the factory. Therefore, we provide a warranty period of 24 months from the date of purchase for this product. **Promed** products conform to their description and the relevant specifications; it is your responsibility to ensure that the product which you have purchased is designed for your intended use.

- In the event of material defects and/or manufacturing errors which appear during proper use and are detected during the warranty period, we will replace all defective parts of the product free of charge within the warranty period, including the labour costs for warranty repairs.
- **The warranty does not cover:**
  - Normal wear on the product
  - Defects which were caused during transportation or storage of the product
  - Defects or damage which were caused through improper use or inadequate maintenance
  - Damage resulting from disregarding the notes in the operating manual
  - Damage resulting from modifications to the product which were not carried out by **Promed**.
  - Damage caused by sharp objects or as a result of torsion, compression, being dropped, an abnormal impact or other actions which are beyond **Promed**'s reasonable control.
  - Parts which are subject to wear (e.g. moving parts such as ball bearings, etc., fasteners, etc.) are generally excluded from the warranty.
- **A claim under the warranty is not enforceable if:**
  - the product is not sent back in its original packaging or in adequately secure packaging;
  - it has been modified or repaired by a person or company other than **Promed** or one of **Promed**'s authorised dealers;
  - the product has been repaired using parts which are not approved by **Promed**;
  - the serial number / LOT number of the product has been removed, deleted, changed or made illegible.
- For reasons of hygiene, any products which have been exposed to direct contact with a body or with bodily fluids (e.g. blood) must be packaged in an additional plastic bag before they are sent back. In this case, a specific note of this situation must be made in the package or in the accompanying documents.
- For components or products which are replaced during the warranty period, a warranty is only provided for the remainder of the original warranty period; provided that the product or part replacement was carried out by **Promed** or one of **Promed**'s authorised dealers.
- The warranty period begins on the day of purchase. Claims under the warranty must be established within the warranty period. Complaints arising after the expiry of the warranty period cannot be considered.
  - Under these warranty conditions, the warranty shall only apply if the date of purchase can be proven by means of a sales receipt or similar.
  - We reserve the right to make technical and visual changes as well as changes to the features.
- This warranty is only valid and enforceable in the country in which it was purchased by the original buyer, provided that it was **Promed**'s intention that the product should be offered for sale in this country. This warranty is also enforceable in any country within the European Economic Area in which **Promed** has an authorised importer or sales partner. Depending on each country, specific and differing warranties and guarantees may be applicable owing to the particular applicable legislation. These legal provisions are neither excluded nor limited by these warranty conditions. To the extent permissible under national law, the warranty period shall not be extended, renewed or otherwise affected by any subsequent resale, repair or replacement of the product.
- The provisions of the UN Convention on Contracts for the International Sale of Goods shall not apply.
- Our warranty conditions do not affect the seller's legal warranty obligation.
- To the fullest extent possible under the applicable mandatory legislation, these warranty conditions constitute your sole and exclusive legal remedy and shall apply in place of any other express or implied warranty conditions. **Promed** is not liable for any unusual, incidental, punitive or consequential damages including, but not limited to loss of profits, loss of use, loss of revenue, costs of replacement equipment or facilities, third party insurance claims, and damage to property which are the result of the purchase or use of the product or occur owing to a breach of warranty, a breach of contract, negligence, product defects or any other legal or statutory circumstances, even though **Promed** was aware of the possibility of such damages. **Promed** is not liable for any delay in warranty claims.
- **Promed** cannot be held liable for any possible translation errors.
- The following information is vital for problem-free processing:
  1. Original sales slip/receipt or dealer's stamp with the date of purchase
  2. Defect detected
  3. Unit name / type / serial / batch number

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**EN** | The requirements of the EC Directive 93/42/EEC in conjunction with EC Directive 2007/47/EC for class IIa medical products has been fulfilled according to the Standard EN 60601-1-2.

CE according to the EC Directive

2014/30/EU and EC Directives 2011/65/EU, 1907/2006/EC, 2005/69/EC and 93/68/EEC.

The CE marking for the appliance refers to the EC Directive 93/42/EEC in conjunction with EC Directive 2007/47/EC.

Appliance classification: Class IIa.

Electromagnetic disturbances are neither emitted from the appliance nor affect the appliance.

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Only transport in the original packaging.

#### Transportation/storage conditions

Temperature:  
-10 °C – +55 °C  
(+14 °F to +131 °F)

Relative humidity:  
10 % - 90 % non-condensing  
Atmospheric pressure:  
700 to 1 060 hPa

#### Operating conditions:

Temperature:  
+5 °C – +40 °C  
(+41 °F to +104 °F)

Relative humidity:  
30 % to 75 % non-condensing  
Atmospheric pressure:  
700 to 1 060 hPa

Service life of the device: 3 years  
Service life of battery: new batteries will last for approx. one week (when used for once a day, TENS, P06 50 mA intensity)  
Service life of electrode pads:  
Pads can be cleaned and reused for up to 10-15 times.

An abrupt change in temperature is to be avoided as condensation may arise in the unit. Only operate the unit after the temperature has been equalized.

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SE

FI

- This device should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this device should be observed to verify normal operation in the configuration in which it will be used.
- Use of accessories other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- When the operating environment is relatively dry, strong electromagnetic interference usually occurs. At this time, the device may be affected as follows:
  - the device stops output;
  - the device turns off;
  - the device restarts;

The above phenomenon does not affect the basic safety and essential performance of the device, and the user can use it according to the instruction. If you want to avoid the above phenomenon, please use it according to the environment specified in the manual.

Table 1

| Declaration - electromagnetic emission  |                |  |
|---|----------------|--|
| The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment. |                |  |
| Emissions test  | Compliance     | Electromagnetic environment - guidance   |
| RF emissions<br>CISPR 11  | Group 1        | The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.              |
| RF emissions<br>CISPR 11  | Class B        | The device is suitable for domestic establishment and in establishment directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions<br>IEC 61000-3-2   | Not applicable |  |
| Voltage fluctuations/<br>flicker emissions<br>IEC 61000-3-3   | Not applicable |  |

Table 2

| Declaration - electromagnetic immunity  |   |   |  |
|---|---|---|--|
| The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment. |   |   |  |
| Immunity test   | IEC 60601 test level  | Compliance level                                    | Electromagnetic environment - guidance   |
| Electrostatic discharge (ESD)<br>IEC 61000-4-2  | ±8 kV contact<br>±2 kV, ±4 kV, ±8 kV,<br>±15 kV air   | ±8 kV contact<br>±2 kV, ±4 kV, ±8<br>kV, ±15 kV air | Floors should be wood, concrete or ceramic tie. If floors are covered with synthetic material, the relative humidity should be at least 30 %.  |
| Electrical fast transient/burst<br>IEC 61000-4-4  | ± 2 kV for power supply lines<br>± 1 kV for input/output lines  | Not applicable                                      | Mains power quality should be that of a typical commercial or hospital environment.  |
| Surge<br>IEC 61000-4-5  | ± 0.5kV, ± 1 kV line(s) to lines<br>± 0.5kV, ± 1 kV, ± 2 kV line(s) to earth  | Not applicable                                      | Mains power quality should be that of a typical commercial or hospital environment.  |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11  | 0 % UT; 0.5 cycle<br>At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°<br><br>0 % UT; 1 cycle and 70 % UT;<br>25/30 cycles<br>Single phase: at 0°<br><br>0 % UT;<br>250/300 cycles | Not applicable                                      | Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000-4-8  | 30 A/m  | 30 A/m  | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.  |
| NOTE: UT is the a.c. mains voltage prior to application of the test level.  |   |   |  |

Table 3

| Declaration - electromagnetic immunity  |   |                  |   |
|---|---|------------------|---|
| The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.   |   |                  |   |
| Immunity test   | IEC 60601 test level  | Compliance level | Electromagnetic environment - guidance  |
| Conducted RF<br>IEC 61000-4-6   | 3V<br>0.15 MHz to 80MHz<br>6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz | Not applicable   | <p>Portable and mobile RF communications equipment should be used no closer to any part of device, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance<br/> <math>d=2.1\sqrt{P}</math> 150 KHz to 80 MHz<br/> <math>d=2.1\sqrt{P}</math> 80 MHz to 800 MHz<br/> <math>d=2.3\sqrt{P}</math> 80 MHz to 2.7 GHz<br/>                     where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
| Radiated RF<br>IEC 61000-4-3  | 10V/m<br>80 MHz to 2.7 GHz  | 10V/m            |   |
| NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.<br>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. |   |                  |   |

- Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which device is used exceeds the applicable RF compliance level above, device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating device.
- Over the frequency range 0.15 MHz to 80 MHz, field strengths should be less than 3 V/m.

Table 4

| Recommended separation distances between portable and mobile RF communications equipment and device  |  |                                      |                                      |
|--|--|--------------------------------------|--------------------------------------|
| The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and device, as recommended below, according to the maximum output power of the communications equipment.  |  |                                      |                                      |
| Rated maximum output power of transmitter<br>W   | Separation distance according to frequency of transmitter<br>m |                                      |                                      |
|  | 0.15 MHz to 80 MHz<br>$d=2.1\sqrt{P}$                          | 80 MHz to 800 MHz<br>$d=2.1\sqrt{P}$ | 80 MHz to 2.7 GHz<br>$d=2.3\sqrt{P}$ |
| 0.01   | 0.12   | 0.12                                 | 0.23                                 |
| 0.1  | 0.38   | 0.38                                 | 0.73                                 |
| 1  | 1.2  | 1.2                                  | 2.3                                  |
| 10   | 3.8  | 3.8                                  | 7.3                                  |
| 100  | 12   | 12                                   | 23                                   |
| For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer<br>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.<br>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. |  |                                      |                                      |

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Shenzhen Dongdixin Technology Co., Ltd.
Shanghai International Holding Corp. GmbH (Europe)
MedNet SWISS GmbH, Bäderstrasse 18, 5400 Baden,